



Compare os 18 modelos | Linha eLight

MODELOS (Potência Total/Potência por Emissor)	Led Azul 450 ± 10 nm	Led Âmbar 590 ± 10 nm	Led Violeta 410 ± 10 nm	Led Vermelho 630 ± 15 nm	Led Infravermelho 850 ± 10 nm	Laser Vermelho 660 ± 10 nm	Laser Infravermelho 808 ± 10 nm	Emissão simultânea
	B	4.160 mW 8 x 520 mW	—	—	—	—	—	—
A	—	800 mW 8 x 100 mW	—	—	—	—	—	—
V	—	—	3.840 mW 8 x 480 mW	—	—	—	—	—
R	—	—	—	1.600 mW 8 x 200 mW	—	—	—	—
I	—	—	—	—	2.080 mW 8 x 260 mW	—	—	—
IV	—	—	1.920 mW 4 x 480 mW	—	1.040 mW 4 x 260 mW	—	—	Sim
RV	—	—	1.920 mW 4 x 480 mW	800 mW 4 x 200 mW	—	—	—	Sim
ARV	—	400 mW 4 x 100 mW	1.920 mW 4 x 480 mW	800 mW 4 x 200 mW	—	—	—	Não
ABR	2.080 mW 4 x 520 mW	400 mW 4 x 100 mW	—	800 mW 4 x 200 mW	—	—	—	Não
AIL	—	800 mW 8 x 100 mW	—	—	—	—	200 mW 2 x 100 mW	Sim
ABV	2.080 mW 4 x 520 mW	400 mW 4 x 100 mW	1.920 mW 4 x 480 mW	—	—	—	—	Não
BR	2.080 mW 4 x 520 mW	—	—	800 mW 4 x 200 mW	—	—	—	Sim
BRL	4.160 mW 8 x 520 mW	—	—	—	—	200 mW 2 x 100 mW	—	Sim
ILB	4.160 mW 8 x 520 mW	—	—	—	—	—	200 mW 2 x 100 mW	Sim
ILV	—	—	3.840 mW 8 x 480 mW	—	—	—	200 mW 2 x 100 mW	Sim
RL	—	—	—	—	—	800 mW 8 x 100 mW	—	—
IL	—	—	—	—	—	—	960 mW 8 x 120 mW	—
IRL	—	—	—	—	—	400 mW 4 x 100 mW	400 mW 4 x 100 mW	Sim